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## 2. Assembly of the spectroscope

The Student Spectroscope in the well known design of Kirchhoff and Bunsen is used for the qualitative observing of spectra. It is equipped with a flint glass prism with a refractive angle of 60°, dispersion  $n_F - n_C = 0.017$ .

An adjustable platform is fitted to the top of a stand, with the equilateral flint glass prism, firmly clamped to it. The following components are also fitted to the platform:

1. The collimator (5) with variable slit (6).

2. The observation tube

is mounted on a swivel holder that can be adjusted by a screw (7). It can be focused for infinity by setting the eyepiece (10).

3. The scale tube (1)

a graduation which is adjusted so that lies in the focal plane of the convergent lens.

The prism is covered by a detachable cap which keeps away interfering light.

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
Determine the lines of a known spectrum (e.g. with a spectrum chart) and mark off the associated wavelengths on the y-axis of a system of coordinates. Mark off the scale reading of the same lines on the x-axis. The connection of these points delivers the diagram showing the interrelationship between scale reading and wavelength.

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For precise measurements it is necessary to use several different spectrums for the calibration in order to get the distances between the single points little.

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1. Keep it away from dust. Always cover the spectroscope with the vinyl dust cover. If dust gets on lenses, blow off with rubber syringe or use camel's hair brush. Optical lenses are very soft, and can easily get scratched by dust particles.

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2. Do not touch the lens. are very difficult to remove. However, if you do, use washed out soft linen moistened with xylol to remove it.

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3. Keep it away from moisture. Moisture is the worst enemy to the optical lenses. Once you get fungus on the lenses, it cannot be returned to the original condition. Keep the instrument in a dry place.

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### **1. The field of view is dark**

The luminous source is dark or dislocated  
 Switch on the source or locate it in front of the slit.

The slit is not opened  
 Open slit

The telescope is swivelled in such an extreme position, that it shows only the invisible spectral region

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Turn the telescope to the left or right.

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The Prism is dislocated  
Remove the cap (4). Loosen the fixing screw (7). Turn the prism as far as the spectrum becomes visible.

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**2. The scale is invisible, obscure  
inclined**

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The scale illumination is not connected or too dull  
Switch on the illumination or change the bulb.

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The slit is dusty  
Clean by means of a sheet of wet paper or a leather cloth after opening the slit completely. Never touch the slit with metal.

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The slit is damaged  
Change the damaged slit